

SPEC SHEET + MEASUREMENTS

SPORT-TEK®

Sport-Tek® Women's Micropique Sport-Wick® Polo. LST650

Product photo



Smooth micropique polos that wick moisture and resist snags.

- 3.8-ounce, 100% polyester tricot
- Snag resistant
- Moisture-wicking
- Double-needle stitching throughout
- Tag-free label
- Taped neck
- Self-fabric collar
- Open placket
- Armhole accent
- Curved back waist seam for flattering fit
- Side vents
- Rental-friendly

CARE INSTRUCTIONS

Machine wash cold with like colors. Only non-chlorine bleach, when needed. Tumble dry low. Remove promptly. Cool iron, if necessary. Do not dry clean.

Sketches



front



back

SPEC SHEET + MEASUREMENTS



Sport-Tek® Women's Micropique Sport-Wick® Polo. LST650

PRODUCT MEASUREMENTS

	XS	S	M	L	XL	XXL	3XL	4XL
Bust	18 1/2	19 1/2	20 1/2	22	23 1/2	25	27	29
Sleeve Length	14 1/8	14 5/8	15 1/8	15 7/8	16 5/8	17 3/8	18 3/8	19 3/8
Body Length at Back	26	26 1/2	27	28	29	30	30 1/2	31

Bust: Measured across the chest one inch below armhole when laid flat.

Sleeve Length: Start at center of neck and measure down shoulder, down sleeve to hem.

Body Length at Back: Measured from high point shoulder to finished hem at back.

SIZE CHARTS

	XS	S	M	L	XL	XXL	3XL	4XL
Size	2	4/6	8/10	12/14	16/18	20/22	24/26	28/30
Bust	32-34	35-36	37-38	39-41	42-44	45-47	48-51	52-55

SPEC SHEET + MEASUREMENTS



Sport-Tek® Women's Micropique Sport-Wick® Polo. LST650

HOW TO MEASURE

BUST



Measure under the arm and around the fullest part of the bust with arms down, keeping tape horizontal.

COLORS

Please note: PMS code information is shown only when an exact match is available.



Black
PMS BLACK
C



Blue Lake
PMS 2129C



Deep Orange
PMS 173C



Deep Red
PMS 201C



Forest Green
PMS 553C



Gold
PMS 1235C



Grey Concrete
PMS COOL
GRAY 9C



Iron Grey
PMS 7540C



Kelly Green
PMS 341C



Light Pink
PMS 3568C



Lime Shock
PMS 375C



Maroon
PMS 504C



Pink
Raspberry
PMS 7425C



Purple
PMS 7672C



Tropic Blue
PMS 7474C



True Navy
PMS 533C



True Red
PMS 200C



True Royal
PMS 7687C



White